

# Exploring CD4+ T Helper Cell Differentiation in Tumor and Peripheral Blood of Clear Cell Renal Carcinoma Patient Using a 50-Plus-Parameter CyTOF Panel

Ling Wang, Shakir Hasan, William Wei, Katrina Thomson, David Howell, Stephen Li, Michael Cohen, Christina Loh  
Standard BioTools Canada Inc., Markham, ON, Canada.

CYTO 2025



## Introduction

CD4+ T helper (Th) cells are a diverse subset of T cells that orchestrate immune responses to acute/chronic infection and tumor growth and are involved in the pathogenesis of autoimmune diseases. Triggered by polarizing cytokines, Th cells differentiate into a range of heterogeneous populations with distinct expression patterns of surface receptors, cytokines and master transcription factors (TFs). Differentiated Th cells such as Th1, Th2, regulatory (Treg), follicular helper (Tfh) and Th17 play discrete roles in combating or promoting specific pathogens or tumors. Recent studies showed that the balance of positive and negative effects of Th cells can influence overall success of cancer immunotherapy. Thus, a deeper understanding of Th cell differentiation and functional states in cancer patients is required to harness the full potential of the immune system to sustain a durable, robust antitumor response.

In this study, we evaluated Th cell subsets in the tumor tissue and peripheral blood from a clear cell renal carcinoma (ccRCC) patient using a comprehensive 50-plus-parameter CyTOF™ panel. This panel includes a diverse array of phenotypic and functional markers identifying Th cell subsets as well as evaluating T cell memory, activation, proliferation, differentiation and exhaustion. Heterogeneous Th cell subsets with divergent composition and functions were observed in both tumor and blood of the ccRCC patient.

Unlike fluorescence-based cytometry, CyTOF technology has low signal spillover and no autofluorescence, and therefore spectral compensation and unmixing are not required. As a result, CyTOF technology enables in-depth single-cell analysis with exceptional resolution and fast panel design. Furthermore, antibody cocktails and stained samples can be frozen for later use and acquisition, enabling a streamlined and flexible workflow in translational and clinical research.

## Materials and methods




## Conclusions

- Deep single-cell profiling using a high-parameter CyTOF panel unveiled an immunosuppressive TME (tumor microenvironment) composed of exhausted T cells, NK cells and suppressive macrophages
- Further in-depth profiling of CD4+ Th cell subsets from a ccRCC patient revealed Th17/Th1 dominant memory phenotypes in the PBMC sample whereas a Th1-polarized exhausted effector phenotype was found in the TDC sample
- Exceptional signal resolution of TFs disclosed a unique correlation between FoxP3 and Bcl-6 in TDC, implying the involvement of Tfr (follicular regulatory T cells) in the TME of ccRCC


## Results

### Deep cellular phenotyping with CyTOF technology reveal bias toward CD4+ T cell infiltration in the tumor microenvironment (TME)





### Comprehensive single-cell profiling unveils an immunosuppressive TME composed of activated and exhausted T cells, NK cells and suppressive macrophages

#### CD45+ cells from ccRCC patient



### In-depth profiling of CD4+ T helper cells shows dominant Th1 polarization in the TME with activated and exhausted effector phenotypes

#### CD4+T cells of ccRCC patient



### Diverse and interrelated cytokine profiling of Th cells confirms Th17/Th1 polarization in ccRCC PBMC

#### PMA/i stimulated ccRCC PBMC

