

Spatial Profiling of Exhausted T Cells Using High-Plex Imaging Mass Cytometry

Thomas D. Pfister¹, Jyh Yun Chwee¹, Qanber Raza¹, Nikesh Parsotam¹, David Howell², Liang Lim¹, Christina Loh¹

¹Standard BioTools Canada Inc., Markham, ON, Canada ²Standard BioTools Inc., South San Francisco, CA, United States

Introduction

The spatial organization and cellular composition of the tumor microenvironment (TME) has the potential to inform clinical and translational researchers about mechanisms of disease progression and therapeutic success. Of particular interest are immune cells, especially T cells, which can become exhausted due to chronic stimulation. This poses a significant challenge in cancer therapy as exhausted T cells have reduced effector functions and sustained expression of inhibitory receptors such as PD-1 and CTLA-4, resulting in failure to effectively eliminate tumor cells.

Imaging Mass Cytometry[™] (IMC[™]) technology is a high-plex imaging technique that enables deep characterization of the heterogeneity and complexity of the TME. The Hyperion[™] XTi Imaging System utilizes IMC technology to provide signal intensities over a wide dynamic range and entails one-step detection of 40-plus markers without issues of tissue autofluorescence, making it ideally suited for spatial biology applications. Whole slide imaging modes and an automated slide loader function enable a streamlined, versatile, scalable workflow for high-throughput analysis.

Methods and materials

The 41-marker panel used in this study was created by adding commercially available expansion panels and single antibodies to the Human Immuno-Oncology IMC Panel, 31 Antibodies. This expands our ability to conduct comprehensive high-plex tumor and immune cell profiling. Whole tumor tissue sections were stained with this comprehensive antibody panel. Tissue Mode imaging of whole slide tumor sections, combined with pixel-clustering analysis, provided a spatially resolved quantitative assessment of specific tumor and immune components within the TME. This approach was further enhanced by a quick tissue scan using Preview Mode, which was used to guide single-cell analysis of selected regions of interest (ROIs) in serial tissue sections that were acquired at single-cell resolution using Cell Mode. Together, these methods successfully delivered quantitative spatial biology analyses.



Figure 1. Pathologist-verified 41-marker antibody panel. The Human Immuno-Oncology IMC Panel is designed to explore immunooncological processes in human tumors. It includes 31 pathologist-verified antibodies in the base panel and is optimized for FFPE tissues. The panel's modular structure allows for customization, making it suitable for various translational and clinical samples. When combined with the Human T Cell Exhaustion IMC Panel and the Maxpar[™] IMC Cell Segmentation Kit, it enables the detection of immune cell subtypes, tumor characteristics and microenvironment components, and the presence of cancer-associated fibroblasts (CAFs). This comprehensive approach enhances understanding of the TME in immuno-oncology research.



Figure 2. Imaging Mass Cytometry workflows. (A) IMC technology offers a

streamlined workflow that simplifies translational and clinical application of multiplexed tissue analysis. The five-step process, which consists of obtaining metalconjugated antibodies, staining tissues with antibody cocktails, imaging tissues with the Hyperion XTi Imaging System, and the collection and analysis of high-dimensional data, can be accomplished in as little as 72 hours (two slides with two 4 mm² ROI each). B) The novel WSI modes for IMC technology offer a customized workflow for specific imaging applications. Here we highlight two simple ways for a user to get started. For single-cell analysis, start with Preview Mode, which provides a rapid scan of the whole tissue and highlights all your stained markers. This helps guide ROI placement to capture single cell-resolution image data using Cell Mode. For pixel-clustering analysis of an entire tissue section, users can first identify the placement of tissue using the rapid Brightfield Mode, followed by the novel Tissue Mode, which generates a high-quality scan of the entire tissue section in a matter of hours with higher spot-size ablations enabling entire tissue analysis using pixel-clustering analysis. Combining these new workflows with the newly available slide loader for the Hyperion XTi Imaging System streamlines IMC application and makes it a useful resource for high-throughput clinical and translational studies.

Conclusions

This work characterizing T cell exhaustion markers in multiple cancers showcases the capabilities of IMC technology and establishes it as a reliable high-plex, highthroughput spatial biology imaging platform. IMC technology is ideally suited for developing future translational and clinical applications and has the potential to help guide personalized therapeutic strategies for cancer treatment.

Results: IMC analysis revealed striking heterogeneity with distinct tumor and immune-rich niches in the TME of the cancer tissues.





(B) The tissue is further divided into 12 clusters following unsupervised pixel-clustering analysis along with hierarchical clustering using the MCD SmartViewer analysis pipeline. Three ROIs were investigated further in (C) and (D) using a serial section of the same tissue block (i, ii and iii). (C) The three ROIs identified in (B) were acquired using Cell Mode and visualized using MCD Viewer (i', ii' and iii'). Scale bar is 200 µm.

Patent and License Information: www.standardbio.com/legal/notices. Trademarks: www.standardbio.com/legal/trademarks. Any other trademarks are the sole property of their respective owners. ©2025 Standard BioTools Inc. All rights reserved.



AACR 2025





(B) The tissue is further divided into 10 clusters following unsupervised pixel-clustering analysis along with hierarchical clustering using the MCD SmartViewer analysis pipeline. Three ROIs were investigated further in (C) and (D) using a serial section of the same tissue block (i, ii and iii). (C) The three ROIs identified in (B) were acquired using Cell Mode and visualized using MCD Viewer (i', ii' and iii'). Scale bar is 200 µm.



