



Flexible and Comprehensive Phenotyping of Tumor and Peripheral Blood Mononuclear Cells in Endometrial Carcinoma

Jessica Elmore, Deeqa Mahamed, Geneve Awong, Gary Impey
Standard BioTools Inc., Markham, ON, Canada

AACR 2025 | Poster #3974

Introduction

Cancer biomarkers have revolutionized management and treatment of the disease, leading to remarkable advancements in personalizing medicine and determining optimal therapeutic combinations. However, researchers face significant challenges in identifying biomarkers in the context of extreme biological heterogeneity, a common feature of many cancers. While immune profiling of peripheral blood is a common approach, many critical disease processes are only evident within the tumor environment. These processes can provide valuable prognostic and diagnostic insights or reveal therapeutic targets that are not detectable in the peripheral immune system.

Endometrial carcinoma (EC) is the most prevalent form of uterine cancer, with its incidence rising in developed countries due to factors such as population aging and increasing obesity rates. Despite a relatively favorable prognosis, with an 80% survival rate at 5 years post-diagnosis, the primary curative treatment remains the total removal of the uterus, ovaries and fallopian tubes. This underscores the need for identifying biomarkers that can lead to more personalized medicine approaches and the development of immunotherapies, particularly for aggressive EC subtypes.

Mass cytometry is a high-plex proteomic technology that simultaneously resolves phenotypic and functional markers, enabling researchers to implement large-scale immunophenotyping strategies that span biological heterogeneity and sample types. A combination of immune and functional profiling is key to elucidating disease mechanisms and revealing predictive biomarkers. Mass cytometry conducted on the CyTOF™ XT PRO system uniquely enables higher-parameter, precise immunophenotyping at greater throughput without the data artifacts introduced by compensation and spectral deconvolution. The ability to easily and rapidly design and modify 50-plus-marker panels, along with flexible sample staining and acquisition workflows and the use of sample multiplexing, makes mass cytometry the premier choice for large and complex clinical studies and drug discovery programs.

Objective

To provide a means to achieve deep phenotyping and functional characterization of multiple clinical samples in a single tube by leveraging modular, ready-to-use Flex-Fit™ panels and the enhanced throughput of the CyTOF XT PRO system. By showcasing characterization of both immune and non-immune cells in tumor tissue and PBMC, this method highlights a rapid, high-plex workflow that provides valuable insights for cancer research and potential therapeutic targets from minimal sample amounts.

Key advantages of mass cytometry on the CyTOF XT PRO system for clinical research

- Capturing phenotypic and functional variation in a single CyTOF panel generates unique biomarkers that can reveal mechanisms of disease activity, drug response and prognostic potential
- **Pre-optimized modular panels** combined with the enhanced throughput of the CyTOF XT PRO system enables fast implementation of large-scale immune profiling studies
- **Sample multiplexing** using a variety of barcoding reagents (Pd, CD45, TeMaI) is a powerful method to harmonize sample sets, reduce batch effects, and improve standardization in multi-site and longitudinal studies
- The CyTOF XT PRO system addresses regulatory requirements with **21 CFR Part 11 compliance-enabling software** ensuring user management, user audit trails and integrity of output files

Materials and methods

Host	Target	Clone	Host	Staining Protocol	Part Number
Human T8Hx40 CyTOF Panel, 3 Antibodies (201308)	CD3	UCHL1	1D6	Surface	2148001
	CD8	SK1	1D6	Surface	2148002
	CD19	SK6	1D6	Surface	2148003
	CD4	SK1	1D6	Surface	2148004
	CD80	SK1	1D6	Surface	2148005
	CD45	SK1	1D6	Surface	2148006
	CD137	SK1	1D6	Surface	2148007
	CD138	SK1	1D6	Surface	2148008
	CD139	SK1	1D6	Surface	2148009
	CD140	SK1	1D6	Surface	2148010
Human T Cell Proliferation CyTOF Panel, 10 Antibodies (201348)	CD3	UCHL1	1D6	Surface	2148011
	CD8	SK1	1D6	Surface	2148012
	CD19	SK6	1D6	Surface	2148013
	CD4	SK1	1D6	Surface	2148014
	CD80	SK1	1D6	Surface	2148015
	CD45	SK1	1D6	Surface	2148016
	CD137	SK1	1D6	Surface	2148017
	CD138	SK1	1D6	Surface	2148018
	CD139	SK1	1D6	Surface	2148019
	CD140	SK1	1D6	Surface	2148020
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148021
	PD-1	29F.1A12	1D6	Surface	2148022
	LAG-3	17A12	1D6	Surface	2148023
	CD3	UCHL1	1D6	Surface	2148024
	CD8	SK1	1D6	Surface	2148025
	CD19	SK6	1D6	Surface	2148026
	CD4	SK1	1D6	Surface	2148027
	CD80	SK1	1D6	Surface	2148028
	CD45	SK1	1D6	Surface	2148029
	CD137	SK1	1D6	Surface	2148030
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148031
	PD-1	29F.1A12	1D6	Surface	2148032
	LAG-3	17A12	1D6	Surface	2148033
	CD3	UCHL1	1D6	Surface	2148034
	CD8	SK1	1D6	Surface	2148035
	CD19	SK6	1D6	Surface	2148036
	CD4	SK1	1D6	Surface	2148037
	CD80	SK1	1D6	Surface	2148038
	CD45	SK1	1D6	Surface	2148039
	CD137	SK1	1D6	Surface	2148040
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148041
	PD-1	29F.1A12	1D6	Surface	2148042
	LAG-3	17A12	1D6	Surface	2148043
	CD3	UCHL1	1D6	Surface	2148044
	CD8	SK1	1D6	Surface	2148045
	CD19	SK6	1D6	Surface	2148046
	CD4	SK1	1D6	Surface	2148047
	CD80	SK1	1D6	Surface	2148048
	CD45	SK1	1D6	Surface	2148049
	CD137	SK1	1D6	Surface	2148050
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148051
	PD-1	29F.1A12	1D6	Surface	2148052
	LAG-3	17A12	1D6	Surface	2148053
	CD3	UCHL1	1D6	Surface	2148054
	CD8	SK1	1D6	Surface	2148055
	CD19	SK6	1D6	Surface	2148056
	CD4	SK1	1D6	Surface	2148057
	CD80	SK1	1D6	Surface	2148058
	CD45	SK1	1D6	Surface	2148059
	CD137	SK1	1D6	Surface	2148060
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148061
	PD-1	29F.1A12	1D6	Surface	2148062
	LAG-3	17A12	1D6	Surface	2148063
	CD3	UCHL1	1D6	Surface	2148064
	CD8	SK1	1D6	Surface	2148065
	CD19	SK6	1D6	Surface	2148066
	CD4	SK1	1D6	Surface	2148067
	CD80	SK1	1D6	Surface	2148068
	CD45	SK1	1D6	Surface	2148069
	CD137	SK1	1D6	Surface	2148070
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148071
	PD-1	29F.1A12	1D6	Surface	2148072
	LAG-3	17A12	1D6	Surface	2148073
	CD3	UCHL1	1D6	Surface	2148074
	CD8	SK1	1D6	Surface	2148075
	CD19	SK6	1D6	Surface	2148076
	CD4	SK1	1D6	Surface	2148077
	CD80	SK1	1D6	Surface	2148078
	CD45	SK1	1D6	Surface	2148079
	CD137	SK1	1D6	Surface	2148080
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148081
	PD-1	29F.1A12	1D6	Surface	2148082
	LAG-3	17A12	1D6	Surface	2148083
	CD3	UCHL1	1D6	Surface	2148084
	CD8	SK1	1D6	Surface	2148085
	CD19	SK6	1D6	Surface	2148086
	CD4	SK1	1D6	Surface	2148087
	CD80	SK1	1D6	Surface	2148088
	CD45	SK1	1D6	Surface	2148089
	CD137	SK1	1D6	Surface	2148090
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148091
	PD-1	29F.1A12	1D6	Surface	2148092
	LAG-3	17A12	1D6	Surface	2148093
	CD3	UCHL1	1D6	Surface	2148094
	CD8	SK1	1D6	Surface	2148095
	CD19	SK6	1D6	Surface	2148096
	CD4	SK1	1D6	Surface	2148097
	CD80	SK1	1D6	Surface	2148098
	CD45	SK1	1D6	Surface	2148099
	CD137	SK1	1D6	Surface	2148100
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148101
	PD-1	29F.1A12	1D6	Surface	2148102
	LAG-3	17A12	1D6	Surface	2148103
	CD3	UCHL1	1D6	Surface	2148104
	CD8	SK1	1D6	Surface	2148105
	CD19	SK6	1D6	Surface	2148106
	CD4	SK1	1D6	Surface	2148107
	CD80	SK1	1D6	Surface	2148108
	CD45	SK1	1D6	Surface	2148109
	CD137	SK1	1D6	Surface	2148110
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148111
	PD-1	29F.1A12	1D6	Surface	2148112
	LAG-3	17A12	1D6	Surface	2148113
	CD3	UCHL1	1D6	Surface	2148114
	CD8	SK1	1D6	Surface	2148115
	CD19	SK6	1D6	Surface	2148116
	CD4	SK1	1D6	Surface	2148117
	CD80	SK1	1D6	Surface	2148118
	CD45	SK1	1D6	Surface	2148119
	CD137	SK1	1D6	Surface	2148120
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148121
	PD-1	29F.1A12	1D6	Surface	2148122
	LAG-3	17A12	1D6	Surface	2148123
	CD3	UCHL1	1D6	Surface	2148124
	CD8	SK1	1D6	Surface	2148125
	CD19	SK6	1D6	Surface	2148126
	CD4	SK1	1D6	Surface	2148127
	CD80	SK1	1D6	Surface	2148128
	CD45	SK1	1D6	Surface	2148129
	CD137	SK1	1D6	Surface	2148130
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148131
	PD-1	29F.1A12	1D6	Surface	2148132
	LAG-3	17A12	1D6	Surface	2148133
	CD3	UCHL1	1D6	Surface	2148134
	CD8	SK1	1D6	Surface	2148135
	CD19	SK6	1D6	Surface	2148136
	CD4	SK1	1D6	Surface	2148137
	CD80	SK1	1D6	Surface	2148138
	CD45	SK1	1D6	Surface	2148139
	CD137	SK1	1D6	Surface	2148140
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148141
	PD-1	29F.1A12	1D6	Surface	2148142
	LAG-3	17A12	1D6	Surface	2148143
	CD3	UCHL1	1D6	Surface	2148144
	CD8	SK1	1D6	Surface	2148145
	CD19	SK6	1D6	Surface	2148146
	CD4	SK1	1D6	Surface	2148147
	CD80	SK1	1D6	Surface	2148148
	CD45	SK1	1D6	Surface	2148149
	CD137	SK1	1D6	Surface	2148150
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148151
	PD-1	29F.1A12	1D6	Surface	2148152
	LAG-3	17A12	1D6	Surface	2148153
	CD3	UCHL1	1D6	Surface	2148154
	CD8	SK1	1D6	Surface	2148155
	CD19	SK6	1D6	Surface	2148156
	CD4	SK1	1D6	Surface	2148157
	CD80	SK1	1D6	Surface	2148158
	CD45	SK1	1D6	Surface	2148159
	CD137	SK1	1D6	Surface	2148160
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148161
	PD-1	29F.1A12	1D6	Surface	2148162
	LAG-3	17A12	1D6	Surface	2148163
	CD3	UCHL1	1D6	Surface	2148164
	CD8	SK1	1D6	Surface	2148165
	CD19	SK6	1D6	Surface	2148166
	CD4	SK1	1D6	Surface	2148167
	CD80	SK1	1D6	Surface	2148168
	CD45	SK1	1D6	Surface	2148169
	CD137	SK1	1D6	Surface	2148170
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148171
	PD-1	29F.1A12	1D6	Surface	2148172
	LAG-3	17A12	1D6	Surface	2148173
	CD3	UCHL1	1D6	Surface	2148174
	CD8	SK1	1D6	Surface	2148175
	CD19	SK6	1D6	Surface	2148176
	CD4	SK1	1D6	Surface	2148177
	CD80	SK1	1D6	Surface	2148178
	CD45	SK1	1D6	Surface	2148179
	CD137	SK1	1D6	Surface	2148180
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148181
	PD-1	29F.1A12	1D6	Surface	2148182
	LAG-3	17A12	1D6	Surface	2148183
	CD3	UCHL1	1D6	Surface	2148184
	CD8	SK1	1D6	Surface	2148185
	CD19	SK6	1D6	Surface	2148186
	CD4	SK1	1D6	Surface	2148187
	CD80	SK1	1D6	Surface	2148188
	CD45	SK1	1D6	Surface	2148189
	CD137	SK1	1D6	Surface	2148190
Human Immune Checkpoint Core CyTOF Panel, 3 Antibodies (201345)	CTLA-4	9A12	1D6	Surface	2148191
	PD-1	29F.1A12	1D6	Surface	2148192
	LAG-3	17A12	1D6	Surface	2148193
	CD3	UCHL1	1D6	Surface	2148194
	CD8	SK1	1D6	Surface	2148195
	CD				