

Targeted Next–Generation Sequencing for Canine Genotyping using the Fluidigm® D3™ Assay Design Engine with the Access Array™ System

Min Lin^{1*}, Mark Neff², Alison Ruhe³, Jiashun Zheng⁴, Nico Tuason¹, Michael Lin¹, Janet Ziegler¹, Honrado Lopez¹ and Julian Walker¹

¹Fluidigm Corporation, 7000 Shoreline Court, Suite 100, South San Francisco, CA 94080
²Van Andel Research Institute, 333 Bostwick Ave NE, Grand Rapids, MI 49503
³ProjectDog, 636 San Pablo Ave, Albany, CA 94706
⁴University of California, San Francisco, CA 94080

Introduction

For proof-of-concept, Fluidigm® D3™ assay design engine and Access Array™ system were implemented for targeted next-generation sequencing for non-human genotyping identification to achieve a rapid and more economical workflow. We assayed 48 variants mined from the genome of a dog affected with adult onset deafness. A prior genome-wide association study (GWAS) had mapped the causative gene to a ~4 Mb interval. The top 48 single-nucleotide polymorphisms (SNPs) from within this interval were converted to assays suitable for NGS analysis. The Fluidigm® D3™ assay design engine was applied to generate target-specific primers, and the Fluidigm® Access Array™ System provided a fast and automated target-enrichment sample preparation. All 48 variants within the 48 samples were successfully genotyped in a single sequencing run. 24 heterozygous variants were excluded because the inherited deafness is a Mendelian recessive defect. Among the remaining 24 variants, three showed low “carrier” rates among a population sampling, consistent with a recent and private mutation in the breed isolate. Our results demonstrate this approach is faster (<three weeks) and more cost effective (<\$0.30/dog/variant) than current platforms for SNP genotyping by NGS.

Study Background




- Inherited deafness:** An adult onset deafness segregates in the Rhodesian Ridgeback breed with an apparent autosomal recessive mode. Affected dogs completely lose the ability to hear between six months and one year of age. Male dogs appear to lose their hearing earlier than female dogs, though there does not appear to be a sex bias of affected status.
- Prior GWAS:** A previous case-control association study (26:230) mapped the cause of the defect to a single locus spanning 4 Mb. There were no canine orthologs of previously-implicated human deafness genes found in the mapped interval.
- Prior whole genome sequencing:** The whole genome of an affected dog was sequenced at 32X coverage using multiple lanes on an Illumina HiSeq. Conventional NGS pipelines for alignment and variant calling on the Amazon Elastic Compute Cloud identified a total of 10,351 sequence variants within the mapped interval.

Conclusions


- This study demonstrates the utility of a simplified, targeted resequencing method for custom SNP genotyping. The results of this experiment, with 48 SNPs typed on 48 dogs, excluded half the variants as being causal and prioritized three SNPs as being strongest candidates.
- 95.0% of NGS reads aligned to the canine reference, and 99.9% of those mapped to the targeted amplicon sites, demonstrating the high quality of the NGS library.
- One of the advantages of genotyping by using the Fluidigm® Access Array™ System is SNP discovery. In this study, in addition to the targeted SNPs, we identified 56 additional SNPs in our 48 target regions.
- The Fluidigm® systems leveraged here provided rapid and cost-effective NGS library preparation, agnostic to the sample type and scale of genetic study by NGS. High read coverage per variant assayed suggests that a greater depth of pooling may be obtained in the future, further increasing the efficiency of this approach and decreasing the cost.

Results


Fluidigm® Systems and Workflow




Primer Validation and Sample QC



Data Analysis Workflow



Data Analysis Workflow



Marker Selection

30 variants were selected for possible causality

10,351 total variants discovered in the affected case genome

4,928 variants homozygous, as expected for a Mendelian recessive defect

79 variants unique to the Rhodesian Ridgeback breed, compared to 30 other dog genomes analyzed (no other breeds are known to exhibit a simple adult onset deafness)

30 variants tied to non-zero phastcons scores in the reference genome sequence, as expected for base pair coordinates of functional effect

0 variants deemed to be non-synonymous coding changes

Analysis of Genotype Data

